

# UX SERIES

UX10, UX12, & UX15 Patient Monitors



High-Resolution Gesture Touch Screen



Bi-directional Communications with Central Monitoring System



HDMI Output for Larger Screen Viewing

# Introducing The UX Series

The Edan UX Series Patient Monitor is a bedside and continuous monitoring solution with flexible mounting solutions for dental care, outpatient settings, and primary care facilities. This advanced patient monitor features a high resolution touchscreen, customizable shortcut keys and a sleek, ergonomic design.


Additionally, the UX Series utilizes iSEAP algorithm for ECG monitoring, which detects up to 33 types of arrhythmias. Blood pressure measurement is quick and accurate with the iFAST algorithm, and the SpO<sub>2</sub> iMAT algorithm offers outstanding motion resistance. With advanced monitoring capabilities, high-capacity storage, and CMS integration, the Edan UX Series is ideal for any environment.

## Features:

- High resolution display & touchscreen
- No-fan design
- 3 hours of battery life (6 hours with upgraded battery)
- 5G WiFi connectivity
- USB, HDMI ports
- Multi-lead analysis
- 3-level alarms with different colors
- Detects up to 33 types of arrhythmia
- Standard 3/5-Lead
- Ample internal storage for trends and event records
- Compatible with MT-207 trolley and wall mount
- iMAT algorithm
- Specialized SEIMP algorithm
- iCUFFS algorithm

**Standard Parameters:** 3/5-lead ECG, HR, RESP, SpO<sub>2</sub>, NIBP, PR, 2-TEMP, IBP

**Optional Features:** Printer, CO<sub>2</sub>, Cardiac Output, F300 Filac Module, Dual IBP



# What's the Difference?

## UX10 Patient Monitor

- ✓ WiFi
  - ✗ Optional CO<sub>2</sub>
  - ✓ 10.1" Touch Screen
  - ✗ Optional Dual IBP Slots
  - ✓ Optional CMS Connectivity
- 

## UX12 Patient Monitor

- ✓ WiFi
  - ✓ Optional CO<sub>2</sub>
  - ✓ 13.3" Touch Screen
  - ✓ Optional Dual IBP Slots
  - ✓ Optional CMS Connectivity
- 

## UX15 Patient Monitor

- ✓ WiFi
- ✓ Optional CO<sub>2</sub>
- ✓ 15.6" Touch Screen
- ✓ Optional Dual IBP Slots
- ✓ Optional CMS Connectivity




## Proprietary Algorithms & Technologies



### G2 CO<sub>2</sub> (Sidestream)

- Superior water trap design for accurate monitoring
- iCARB™ algorithm with intelligent CO<sub>2</sub> pseudo wave identification technology
- Sampling rate as low as 50ml/min



### ECG

- Automatic lead type detection
- Industry leading iSEAP™ algorithm with auto-detection of 33 types of arrhythmias
- SEMIP® algorithm with 208 ECG findings over age/gender diversities



### NIBP

- Dual dust filter design means no blockage inside and provides accurate NIBP readings
- Unique cleaning mode for routine maintenance
- iCUFS™ algorithm with smart deflation technology

### SpO<sub>2</sub>

- iMAT algorithm with motion resistance and low perfusion resistance performance
- Reference reading of Perfusion Index (PI) from 0 to 10 according to perfusion changes
- Simultaneous measurements of SpO<sub>2</sub> and NIBP of the same limb

# Configurations

## UX10

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 1-TEMP, PR, HR and internal memory for previous data review.

## UX10.P

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 1-TEMP, PR, HR, internal memory for previous data review & printer installed.

## UX12

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR and internal memory for previous data review.

## UX12.P

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR, internal memory for previous data review & printer installed.

## UX12-G2

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR and internal OEM Edan Sidestream G2 CO<sub>2</sub> module. G2 CO<sub>2</sub> uses traditional water traps and any male luer lock cannula with generic connection. Internal memory for previous data review.

## UX12-G2.P

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR and internal OEM Edan Sidestream G2 CO<sub>2</sub> module. Internal memory for previous data review and printer installed.

## UX15

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR and internal memory for previous data review.

## UX15.P

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR and internal memory for previous data review & printer installed.

## UX15-G2

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR and internal OEM Edan Sidestream G2 CO<sub>2</sub> module. G2 CO<sub>2</sub> uses traditional water traps and any male luer lock cannula with generic connection. Internal memory for previous data review.

## UX15-G2.P

Standard Configuration with Touch Screen, WiFi, 3/5-lead ECG, RESP, NIBP, Edan SpO<sub>2</sub>, 2-TEMP, PR, HR and internal OEM Edan Sidestream G2 CO<sub>2</sub> module. Internal memory for previous data review and printer installed.

## Included Accessories

### STANDARD ACCESSORIES

- ECG cable, 3-lead, snap, AHA, 3.4m — **01.57.471388**
- SpO<sub>2</sub> Finger Sensor, Adult, 2.5m, reusable - direct connect 7 pin — **02.57.225029**
- NIBP Cuff, Adult, 27cm-35cm, reusable — **Cuff.E9**
- NIBP Tube — **01.59.473007**
- Adult skin temperature probe — **01.15.040225**
- Rechargeable Lithium-Ion Battery (10.8V, 2550mAh) — **01.21.064380**

### G2 ACCESSORIES

- Water Trap — **02.01.210520**
- Disposable CO<sub>2</sub> Sampling line with male luer lock — **4410-10-25**
- Adult Nasal CO<sub>2</sub> sampling cannula — **4000-7-25**

## Optional Accessories

### SPO2 SENSORS

- SpO<sub>2</sub> Finger Sensor, Adult, 2.5m, reusable — **SH1.DB9**
- SpO<sub>2</sub> Wrap Sensor, Neonate, 1m, reusable — **SH3.DB9**
- SpO<sub>2</sub> Silicone Soft-tip Sensor, Adult, 1m, reusable — **SH4.DB9**
- SpO<sub>2</sub> Silicone Soft-tip Sensor, Pediatric, 1m, reusable — **SH5.DB9**
- SpO<sub>2</sub> Ear Clip Sensor, Adult/Pediatric, 1m, reusable — **SH6.DB9**
- SpO<sub>2</sub> 7-pin Extension Cable, 2m — **01.57.471068**
- SpO<sub>2</sub> 7-pin Extension Cable, 4m — **01.57.471789**

### CUFFS

- NIBP Cuff, Infant, 10-15cm, reusable — **Cuff.E5**
- NIBP Cuff, Small Child, 13-17cm, reusable — **Cuff.E6**
- NIBP Cuff, Child, 16-21cm, reusable — **Cuff.E7**
- NIBP Cuff, Small Adult, 20.5-28cm, reusable — **Cuff.E8**
- NIBP Cuff, Adult, 27cm-35cm, reusable — **Cuff.E9**
- NIBP Cuff, Large Adult, 34cm-43cm, reusable — **Cuff.E10**
- NIBP Thigh Cuff, Adult, 42cm-54cm, reusable — **Cuff.E11**

### NIBP TUBING

- NIBP Tube (3m) with connector — **01.59.473007**

# Specifications

## CLASSIFICATION

Anti-electroshock type  
Class I equipment and internal powered equipment  
Anti-electroshock degree: CF  
Ingress Protection: IPX1  
Working system: Continuous operation equipment  
Compliant with Standards: IEC 60601-1; IEC 60601-1-2; EN 60601-1; EN 60601-1-2; IEC 80601-2-49

## PHYSICAL SPECIFICATION

UX10  
Dimensions  
268.5mm(W)×209.7mm(H)×170.5mm(D)  
Weight: 3.5 kg

UX12  
Dimensions  
330.1mm(W)×244.2mm(H)×176.2mm(D)  
Weight: 4.0 kg

UX15  
Dimensions  
383.3mm(W)×268.7mm(H)×182.4 mm(D)  
Weight: 5.0 kg

## ENVIRONMENT SPECIFICATION

Temperature  
Working: 0°C to 40°C (32°F to 104°F)  
Transport and Storage: -20°C to 60°C (-4°F to 140°F)  
Relative Humidity  
Working: 15%RH to 95%RH (non-condensing)  
Transport and Storage: 10%RH to 95%RH (non-condensing)  
Barometric Pressure  
Working: 57 kPa to 107.4 kPa  
Transport and Storage: 16 kPa to 107.4 kPa

## POWER SUPPLY

AC Voltage: 100 V to 240 V  
Input Current: 0.6 A to 0.3 A  
Frequency: 50 Hz/60 Hz  
Over Current Fuse Protection: Support

## BATTERY

Operating Time  
UX10: One battery (2550 mAh) ≥ 4 h  
UX12: One battery (5100 mAh) ≥ 8 h  
UX15: One battery (2550 mAh) ≥ 3 h  
One battery (5100 mAh) ≥ 6 h  
Condition  
At 25±5 °C, with new fully charged battery/batteries, 3-lead ECG cable and SpO<sub>2</sub> sensor connected, NIBP module set at an interval of 15 minutes, and screen brightness set to "1".  
Charge Time  
UX10/UX12: One battery (2550 mAh) ≤ 3.5 h (monitor is off)  
One battery (5100 mAh) ≤ 6.5 h (monitor is off)  
UX15: One battery (2550 mAh) ≤ 5.5 h (monitor is on or standby)  
One battery (5100 mAh) ≤ 11 h (monitor is on or standby)  
Condition  
Environment temperature: 20°C to 30°C ≥ 90% charged  
Charge/Discharge Cycle: 300 times

## DISPLAY

UX10  
Screen: 10.1" color TFT (optional touchscreen)  
Resolution: 1024×600  
Waveforms: 8  
One power LED  
One AC power LED  
One alarm LED  
One battery LED

UX12  
Screen: 13.3" color TFT (optional touchscreen)  
Resolution: 1920×1080

Waveforms: 10  
One power LED  
One AC power LED  
One alarm LED  
One battery LED

UX15  
Screen: 15.6" color TFT (optional touchscreen)  
Resolution: 1920×1080  
Waveforms: 12  
One power LED  
One AC power LED  
One alarm LED  
One battery LED

## RECORDER

Record Width: 48 mm  
Record Paper Width: 50 mm  
Paper Speed: 12.5 mm/s, 25 mm/s, 50 mm/s  
Number of waveform channels: A maximum of 3

## DATA STORAGE

Trend Data  
48 hours @ 1 s  
240 hours @ 1 min  
NIBP Measurement: At least 1600 sets  
Alarm Events: Up to 1800 sets  
Full Disclosure Waveform: 48 hours @ 1 s

## WIFI

IEEE: 802.11a/b/g/n  
Bandwidth: 20 M/40 M  
Frequency Band: 2.4 G/5 G  
Modulation Mode  
802.11a/g/n: OFDM  
802.11b: CCK and DSSS  
Antenna  
Single antenna, 2.4 GHz gain ≤ 3.0  
5 GHz gain ≤ 3.5  
Power voltage: 3.3 V  
Bluetooth: Dual mode BLE4.0+EDR/BR  
Connector: SDIO2.0 for Wi-Fi, UART for BT

## ECG

Lead Mode  
3 Electrodes: I, II, III  
5 Electrodes: I, II, III, aVR, aVL, aVF, V  
6 Electrodes: I, II, III, aVR, aVL, aVF, and leads corresponding to Va Vb.  
10 Electrodes: I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6  
Electrode Standard: AHA, IEC  
Display Sensitivity (Gain Selection)  
1.25 mm/mV (×0.125), 2.5 mm/mV (×0.25), 5 mm/mV (×0.5), 10 mm/mV (×1), 20 mm/mV (×2), 40 mm/mV (×4), AUTO gain

Sweep: 6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s  
Bandwidth (-3dB) (Output amplitude relative to that for a 5 Hz sinusoidal input signal)  
Diagnosis: 0.05 Hz to 150 Hz  
ST: 0.05 Hz to 40 Hz  
Monitor and Monitor (Hi-Fi): 0.5 Hz to 40 Hz  
Enhanced: 2 Hz ~18 Hz  
Customized: High-pass Filter and Low-pass Filter  
CMRR (Common Mode Rejection Ratio)  
Diagnosis: > 95 dB  
ST: > 105 dB  
Monitor and Monitor (Hi-Fi): > 105 dB  
Surgery: > 105 dB  
Enhanced: > 105 dB  
Customized: > 105 dB (Low-pass Filter < 40 Hz) > 95 dB (Low-pass Filter > 40 Hz)

## Hum Filter

In diagnosis, ST, monitor, surgery, enhanced and customized modes: 50 Hz/60 Hz (Hum Filter can be turned on or off manually)  
Differential Input Impedance: > 5 MΩ  
Input Signal Range: ±10 mV PP  
Accuracy of Signal Reproduction  
An error of ≤ ±20 % of the nominal value of the output or ±100 μV, whichever is greater.  
The total error and frequency response comply with IEC 60601-2-27: 2011, Sect. 201.12.1.101.1.

Electrode Offset Potential Tolerance: ±800 mV  
Auxiliary Current (Leads off detection)  
Active electrode: < 100 nA  
Reference electrode: < 900 nA  
Recovery Time After Defibrillation  
< 5 s (measured without electrodes as IEC60601-2-27:2011, Sect. 201.8.5.5.1 requires.)

Leakage Current of Patient: < 10 μA  
Scale Signal: 1 mV PP, accuracy is ±5%  
System Noise: < 30 μVPP  
Multichannel Crosstalk  
≤ 5% of the input signal  
Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.5.

Frequency and Impulse Response  
Frequency response: Input a 5 Hz, 1 mV sine wave signal, and the output signal amplitude remains within the range of 71 % to 110 % at 0.67 Hz and 40 Hz. Input a 1 Hz, 1.5 mV 200 ms triangular wave input signal, and the output shall be within 11.25 mm~15 mm.  
Impulse response:  
Displacement value: ≤ 0.1 mV  
Slope: ≤ 0.3 mV/s following the end of the pulse.  
Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.8.

ESU Protection  
Cut mode: 300 W  
Coagulation mode: 100 W  
Restore time: ≤ 10 s  
Electrosurgical Interference Suppression  
Test according to ANSI/AAMI EC13:2002, Sect. 5.2.9.14. Complied with ANSI/AAMI EC13:2002, Sect. 4.2.9.14.

Sampling Frequency: 1000 Hz  
Sampling Channel Switch Time: < 80 μs  
A/D Precision: 24 Bits (Minimum resolution: 0.077uV/LSB)  
Baseline Reset Time: < 3 s  
Pace Pulse without overshoot  
Pulse Indicator

Pulse is marked if the requirements of IEC 60601-2-27: 2011, Sect. 201.12.1.101.12 are met:  
Amplitude: ±2 mV to ±700 mV  
Width: 0.1 ms to 2.0 ms  
Ascending time: 10 μs to 100 μs

Pulse Rejection  
Pulse is rejected if the requirements of IEC 60601-2-27: 2011, Sect. 201.12.1.101.13 are met:  
Amplitude: ±2 mV to ±700 mV  
Width: 0.1 ms to 2.0 ms  
Ascending time: 10 μs to 100 μs

Pace Pulse Detecting Lead: one among I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6  
Heart Rate  
HR Calculation  
Range

ADU: 15 bpm to 300 bpm  
PED/NEO: 15 bpm to 350 bpm  
Accuracy: ±1% or 1 bpm, whichever is greater  
Resolution: 1 bpm (For display and alarm limit)  
Sensitivity: ≥ 300 μVPP  
QRS Detection Range

The detection range has exceeded the requirement described in the standard:  
Width: 70 ms~120 ms for adult, 40 ms~120 ms for Pediatric/neonate.  
Amplitude: 0.5 mv~5 mv  
In adult mode, these two signals are not responded: 1. when QRS amplitude of 0.15 mV or less is applied; 2. when QRS duration of 10 ms and QRS amplitude of 1 mV or less is applied.  
Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.15.

## PVCs Range

ADU: (0 to 300) PVCs/ min  
PED/NEO: (0 to 350) PVCs/ min  
Resolution: 1 PVCs/min (For display and alarm limit)  
Pause/min  
Range: ADU/PED/NEO: (0 to 30) pauses/min  
Resolution: 1 pause/min (For display and alarm limit)  
ST value

Range: -2.0 mV to +2.0 mV  
Accuracy: -0.8 mV to +0.8 mV: ±0.02 mV or 10%, whichever is greater. Beyond this range: not specified.  
Resolution: 0.01 mV (For display and alarm limit)

# Specifications

<p><b>QT measurement</b> Range: 200 ms ~ 800 ms Resolution: 4 ms (For display) Accuracy: ± 30 ms</p> <p><b>QTc measurement</b> Range: 200ms ~ 800 ms Resolution: 1 ms (For display and alarm limit)</p> <p><b>QTc measurement</b> Range: 200ms ~ 800 ms Resolution: 1 ms (For display and alarm limit)</p> <p><b>ΔQTc measurement</b> Range: -600 ms ~ 600 ms Resolution: 1 ms (For display and alarm limit)</p> <p><b>HR Averaging Method</b> Method 1: Heart rate is computed by excluding the minimum and maximum values from the 12 most recent RR intervals and averaging the residual 10 RR intervals. Method 2: If each of three consecutive RR intervals is greater than 1200 ms, then the four most recent RR intervals are averaged to compute the HR.</p> <p><b>Range of Sinus and SV Rhythm Tachy</b> Adult: RR interval for 5 consecutive QRS complex ≤ 0.5 s. Pediatric/neonatal: RR interval for 5 consecutive QRS complex ≤ 0.375 s.</p> <p><b>Normal</b> Adult: 0.5 s &lt; RR interval for 5 consecutive QRS complex &lt; 1.5 s. Pediatric/neonatal: 0.375 s &lt; RR interval for 5 consecutive QRS complex &lt; 1 s.</p> <p><b>Brady</b> Adult: RR interval for 5 consecutive QRS complex ≥ 1.5 s. Pediatric/neonatal: RR interval for 5 consecutive QRS complex ≥ 1 s.</p> <p><b>Range of Ventricular Rhythm V-Tach</b> 5 consecutive ventricular beats and ventricular HR ≥ 100 bpm.</p> <p><b>Vent Rhythm</b> 5 consecutive ventricular beats, and 20 bpm ≤ ventricular HR &lt; 40 bpm.</p> <p><b>Vent Brady</b> 5 consecutive ventricular beats, and ventricular HR &lt; 20 bpm.</p> <p><b>Acc. Vent Rhythm</b> 5 consecutive ventricular beats, and 40 bpm ≤ ventricular HR &lt; 100 bpm.</p> <p><b>Maximum Start-up Alarm Time for Tachycardia V-Tach (Ventricular Tachycardia) 1 mV 206 bpm</b> Gain 0.5: 10 s Gain 1.0: 10 s Gain 2.0: 10 s</p> <p><b>V-Tach (Ventricular Tachycardia) 2 mV 195 bpm</b> Gain 0.5: 10 s Gain 1.0: 10 s Gain 2.0: 10 s</p> <p><b>Response Time of Heart Rate Meter to Change in HR</b> HR range: 80 bpm to 120 bpm Range : Within 11 s HR range: 80 bpm to 40 bpm Range : Within 11 s</p> <p><b>Tall T-wave Rejection</b> Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.17 minimum recommended 1.2 mV T-Wave amplitude</p> <p><b>Accuracy of Heart Rate Meter and Response to Irregular Rhythm</b> Complied with IEC 60601-2-27: 2011, Sect. 201.7.9.2.9.101 b) 4), the HR value after 20 seconds of stabilization is displayed as follows: Ventricular bigeminy: 80 bpm ± 1 bpm Slow alternating ventricular bigeminy: 60 bpm ± 1 bpm Rapid alternating ventricular bigeminy: 120 bpm ± 1 bpm Bidirectional systoles: 91 bpm ± 1 bpm</p> <p><b>Time to Alarm for Heart Rate alarm conditions</b> Asystole alarm: ≤ 10 s HR low alarm: ≤ 10 s HR high alarm: ≤ 10 s</p>	<p><b>Arrhythmia analysis</b> Asystole V-Fib/V-Tach Couplet Vent Rhythm PVC Bigeminy PVC Trigeminy Tachy R on T PVC Irr Rhythm Brady Missed Beat Pacer not Pacing Vent Brady Pacer not Capture VEB Run PVCs Acc. Vent Rhythm IPVC Non-Sustain VT Multiform PVCs Pause/min High Pause Afib PAC Bigeminy PVCs High Low Voltage(Limb) ExtremeBrady PAC Trigeminy Wide QRS Tachy Sustain VT ExtremeTachy V-Tach</p> <p><b>RESP</b> Method: Impedance between RA-LL, RA-LA Measurement lead Options are lead I and II. The default is lead II. Calculation Type: Manual, Automatic Baseline Impedance Range 200 Ω to 2500 Ω (with ECG cables of 1 KΩ resistance) Measuring Sensitivity Within the baseline impedance range: 0.3 Ω Waveform Bandwidth: 0.2 Hz to 3.3 Hz (-3 dB) Respiration Excitation Waveform Sinusoid, 45.6 kHz (±10%), &lt; 350 μA Gain Selection: ×0.25, ×0.5, ×1, ×2, ×3, ×4, ×5 Sweep: 6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s No Breath Detected Alarm Time Setup 10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s; default value is 20 s. Measuring Range: 0 rpm ~ 200 rpm Resolution: 1 rpm (For display and alarm limit) Accuracy ±1 rpm (0 rpm~120 rpm) ±2 rpm (121 rpm~200 rpm) Update time: 1s</p> <p><b>NIBP</b> Technique: Oscillometry Mode: Manual, Auto, Continuous, Sequence Measuring Interval in AUTO Mode (unit: minute) 1/ 2/ 2.5/ 3/ 4/ 5/ 10/ 15/ 30/ 60/ 90/ 120/ 180/ 240/ 360/ 480 and User Define Continuous: 5 min, interval is 5 s Measuring Parameter: SYS, DIA, MAP, PR Pressure Unit: kPa, mmHg, cmH2O Measuring Range Adult Mode SYS: 25 mmHg to 290 mmHg DIA: 10 mmHg to 250 mmHg MAP: 15 mmHg to 260 mmHg Pediatric Mode SYS: 25 mmHg to 240 mmHg DIA: 10 mmHg to 200 mmHg MAP: 15 mmHg to 215 mmHg Neonatal Mode SYS: 25 mmHg to 140 mmHg DIA: 10 mmHg to 115 mmHg MAP: 15 mmHg to 125 mmHg Alarm Type: SYS, DIA, MAP Cuff Pressure Measuring Range 0 mmHg to 300 mmHg Pressure Resolution 1 mmHg (For display and alarm limit) Maximum Mean Error: ±5 mmHg Maximum Standard Deviation: 8 mmHg Maximum Measuring Period Adult/Pediatric: 120 s Neonate: 90 s Typical Measuring Period iCUPS measurement 20 s to 35 s (depend on HR/motion disturbance) (measured with E9 cuff, default inflation value, PR is set as 80 bpm and systolic pressure within 100~120mmHg) iFAST measurement 15 s (depend on SYS, arm circumference and HR/motion disturbance)(measured with E8 cuff, PR is set as 80 bpm and systolic pressure within 100~120 mmHg)</p>	<p>Dual Independent Channel Overpressure Protection Adult: (297±3) mmHg Pediatric: (245±3) mmHg Neonatal: (147±3) mmHg Pre-inflation Pressure Adult: 80/100/120/140/150/160/180/200/220/240 mmHg Pediatric: 80/100/120/140/150/160/180/200 mmHg Neonatal: 60/70/80/100/120 mmHg Venipuncture pressure Adult Default: 80 mmHg Options 20 mmHg, 30 mmHg, 40 mmHg, 50 mmHg, 60 mmHg, 70 mmHg, 80 mmHg, 90 mmHg, 100 mmHg, 110mmHg, 120 mmHg Pediatric Default: 60 mmHg Options 20 mmHg, 30 mmHg, 40 mmHg, 50 mmHg, 60 mmHg, 70 mmHg, 80 mmHg Neonatal Default: 40 mmHg Options: 20 mmHg, 30 mmHg, 40 mmHg, 50 mmHg</p> <p><b>SPO2</b> Measuring Range: 0% to 100% Resolution: 1% (For display and alarm limit) Data Update Period: 1 s Accuracy Adult /Pediatric ±2% (70% to 100% SpO<sub>2</sub>) Undefined (0% to 69% SpO<sub>2</sub>) Neonate ±3% (70% to 100% SpO<sub>2</sub>) Undefined (0% to 69% SpO<sub>2</sub>) Sensor Red Light: (660±3) nm Infrared Light (905±10) nm Emitted Light Energy: &lt;15 mW PI Measuring Range 0.05% to 20%, invalid PI value is -?-, measurement error is not defined. Resolution 0.01% (0.05%-9.99%) 0.1% (10.0%-20.0%)</p> <p><b>PR</b> Parameter PR (SpO<sub>2</sub>) Measuring range: 20 bpm to 300 bpm Accuracy: ±2 bpm Resolution: 1 bpm (For display and alarm limit) PR (NIBP) Measuring range: 40 bpm to 240 bpm Accuracy: ±3 bpm or 3.5%, whichever is greater Resolution: 1 bpm (For display) PR (IBP) Measuring range: 20 bpm to 300 bpm Accuracy: ±2 bpm or ±2%, whichever is greater (30 bpm to 300 bpm); Undefined (20 bpm to 29 bpm) Resolution: 1 bpm (For display and alarm limit)</p> <p><b>TEMP</b> Technique: Thermal resistance Position: Skin, cavity Measure Parameter T1, T2, TD (the absolute value of T2 minus T1) Channel: 2 Sensor Type: YSI-10K and YSI-2.252K Unit: °C, °F Measuring Range: 0 °C to 50 °C (32 °F to 122 °F) Resolution (For display and alarm limit): 0.1 °C (0.1 °F) Accuracy<sup>1</sup> (1) Sensor accuracy: 25 °C~45 °C: ±0.1; Others: ±0.2 °C (2) Without sensor accuracy: ±0.1 °C Refresh Time: Every 1 s to 2 s Temperature Calibration: At an interval of 5 to 10 minutes Measuring Mode: Direct Mode Transient Response Time: ≤ 30 s</p> <p><b>IBP</b> Technique: Direct invasive measurement Channe: 2 channels IBP Measure Measuring Range: (-50 to +360) mmHg Resolution (For display and alarm limit): 1 mmHg (For display and alarm limit) Accuracy (not including sensor): ±2% or ±1 mmHg, whichever is greater</p>
--	---	--

# Specifications

**PPV**  
 Measuring Range: 0%~50%  
 Pressure Unit: kPa, mmHg, cmH<sub>2</sub>O  
 Pressure sensor  
 Sensitivity: 5 μV/V/mmHg  
 Impedance Range: 300 Ω to 3000 Ω  
 Filter: DC~ 12.5 Hz; DC~ 40 Hz  
 Zero Range: ±200 mmHg  
 Pressure Calibration Range  
 IBP (excluding ICP): 80 mmHg to 300 mmHg  
 ICP: 80 mmHg to 300 mmHg  
 Volume Displacement: 7.4 x 10<sup>4</sup> mm<sup>3</sup>/100mmHg

**CO**  
 Technique: Thermodilution Technique  
 Measure Parameters: C.O., TB, TI  
 Measuring Range  
 C.O.: 0.1 L/min to 20 L/min  
 TB: 23 °C to 43 °C (73.4 °F to 109.4 °F)  
 TI: -1 °C to 27 °C (30.2 °F to 80.6 °F)  
 Resolution (For display and alarm limit)  
 C.O.: 0.1 L/min  
 TB, TI: ±0.1 °C (±0.1 °F)  
 Accuracy  
 C.O.: ±5% or ±0.2 L/min  
 TB: ±0.1 °C (±0.18 °F) (not including sensor)  
 TI: ±0.1 °C (±0.18 °F) (not including sensor)

**CO<sub>2</sub>**  
 Intended Patient: Adult, pediatric, neonatal  
 Measure Parameters: etCO<sub>2</sub>, FICO<sub>2</sub>, AwRR  
 Unit: mmHg, %, kPa  
 Measuring Range  
 etCO<sub>2</sub>: 0 mmHg to 152 mmHg  
 FICO<sub>2</sub>: 0 mmHg to 50 mmHg  
 AwRR: 0 rpm to 150 rpm  
 Resolution (For display and alarm limit)  
 etCO<sub>2</sub>: 1 mmHg (For display and alarm limit)  
 FICO<sub>2</sub>: 1 mmHg (For display and alarm limit)  
 AwRR: 1 rpm (For display and alarm limit)  
 Accuracy  
 etCO<sub>2</sub>  
 Measurement conditions: Ambient temperature: (25 ± 3) °C  
 Barometric pressure: (760 ± 10) mmHg  
 Balance gas: N<sub>2</sub>  
 Sample gas flowrate: 100 ml/min  
 I/E ratio: 1:2  
 AwRR ≤ 80 rpm  
 0 mmHg to 40 mmHg  
 ± 2 mmHg  
 41 mmHg to 70 mmHg  
 ± 5% of reading  
 71 mmHg to 100 mmHg  
 ± 8% of reading  
 101 mmHg to 150 mmHg  
 ± 10% of reading  
 AwRR: ±1 rpm  
 Drift of Measure Accuracy  
 Meets the requirements of the measure accuracy  
 Sample Gas Flowrate  
 50 ml/min, 70 ml/min or 100 ml/min (optional), accuracy: ±15 ml/min  
 Warm-up Time  
 Display reading within 20 s; reach to the designed accuracy within 2 minutes.  
 Rise Time  
 < 400 ms (with 2 m gas sampling tube, sample gas flowrate: 100 ml/min)  
 < 500 ms (with 2 m gas sampling tube, sample gas flowrate: 70 ml/min)  
 < 1000 ms (with 2 m gas sampling tube, sample gas flowrate: 50 ml/min)  
 Response Time  
 < 4 s (with 2 m gas sampling tube, sample gas flowrate: 100 ml/min & 70 ml/min)  
 < 5.5 s (with 2 m gas sampling tube, sample gas flowrate: 50 ml/min)  
 Work Mode: Standby (default), measure  
 O<sub>2</sub> Compensation  
 Range: 0% to 100%  
 Resolution: 1%  
 Default: 16%

N<sub>2</sub>O Compensation  
 Range: 0% to 100%  
 Resolution: 1%  
 Default: 0%

AG Compensation  
 Range: 0% to 20%  
 Resolution: 0.1%  
 Default: 0%

Humidity Compensation Method  
 ATPD (default), BTPS

Barometric Pressure Compensation  
 Automatic (The change of barometric pressure will not add additional errors to the measurement values.)

Zero Calibration: Support  
 Calibration  
 Support (It is recommend to be operated by trained personal.)

Alarm: etCO<sub>2</sub>, FICO<sub>2</sub>, AwRR  
 No Breath Detected Alarm Delay  
 10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s;  
 default value is 20 s.

Data Sample Rate: 100 Hz  
 etCO<sub>2</sub> Change<sup>1</sup>  
 With 2 m gas sampling tube, sample gas flowrate: 100 ml/min  
 AwRR ≤ 80 rpm, meet the accuracy mentioned above;  
 AwRR > 80 rpm, etCO<sub>2</sub> descends 8%;  
 AwRR > 120 rpm, etCO<sub>2</sub> descends 10%  
 With 2 m gas sampling tube, sample gas flowrate: 70 ml/min  
 AwRR ≤ 60 rpm, meet the accuracy mentioned above;  
 AwRR > 60 rpm, etCO<sub>2</sub> descends 8%;  
 AwRR > 90 rpm, etCO<sub>2</sub> descends 10%;  
 AwRR > 120 rpm, etCO<sub>2</sub> descends 15%  
 With 2 m gas sampling tube, sample gas flowrate: 50 ml/min  
 AwRR ≤ 20 rpm, meet the accuracy mentioned above;  
 AwRR > 20 rpm, etCO<sub>2</sub> descends 8%;  
 AwRR > 40 rpm, not applicable;