



EDANUSA, CARDIOLINE, QRS DIAGNOSTICS, MDPRO,
MEDITECH, MIR SPIROMETRY PRODUCTS,
MINDRAY DATASCOPE MONITORING




HD3 Bladder Scanner

12 slice 3D image,
no calibration needed



INPUT

Touchable Screen, convenient for the input of patient information.




INTUITION

Double orthogonal images display after the scanning.




CONVENIENCE

USB connecting and printing function.



DISPLAY

Large-sized display for the ultrasonic images of the bladder.



ASSISTANCE

Crosshair projection and automatic outline drawing of the bladder.



MANAGEMENT

Patient data management software.



HD3 BLADDER SCANNER

CALL FOR MORE INFO
888.602.PROS



EDANUSA, CARDIOLINE, QRS DIAGNOSTICS, MDPRO,
MEDITECH, MIR SPIROMETRY PRODUCTS,
MINDRAY DATASCOPE MONITORING



Display

7" TFT LCD (16.9) 800x480

Probe

3D Mechanical Probe
Frequency: 2.5MHz
B Mode Display Angle: 120°

Accuracy

The accuracy result was calculated according to the scanning instruction and on the Dansk Phantom 616 Phantom: 0 to 999 ml \pm 8-10% \pm 10ml

Configuration:

Pad Scan HD3 Main Unit x1
Pad Scan HD3 Base x 1
Battery Charger x1
Battery x2
PC Management Software x1
Aluminum Protective Carrying Case x1
Dimension of the Carrying Case: 440x230x380 (mm)



Others

Double Power Supply Mode
Lithium Battery: 11.IV 2200mAh/AC 100V~240V
Touchable Screen
USB Storage Supported
Multi-Languages Selection

Scanning Time

12 multi-angles ultrasonic images' capturing time <3 seconds



CALL FOR MORE INFO
888.602.PROS

Address: 4204 Jutland Dr, Ste: B, San Diego, CA 92117
Tel. No.: 888-602-7767
Website: <http://www.medicaldevicepros.com>
Email: info@meddevicepros.com