

# SE-1200 Express SE-1200 Express Basic 12-Channel ECG



The Edan SE-1200 Express Basic (EB) Machine boasts exceptional cost performance with both accurate clinical performance as well as user-friendly design.

## What's The Difference?

### SE-1200 Express Basic


- Short-cut keys only
- 100 EKG Memory

### SE-1200 Express

- Touch screen and short cut keys
- 200 EKG memory



Wireless AP



One-Touch Operation

# User-Friendly Design

## Highlighted Features

### SmartECG Viewer



- Take custom measurements
- Compare previous ECGs
- Increase confidence when diagnosing
- Easily export ECG reports

### Data Management

100 ECG internal memory (can be enlarged by USB flash disk)

Optional SmartECG viewer

### Exceptional Clinical Performance

Weak signal detection with 24 bits A/D converter technology

Anti-noise technology with comprehensive filters

High CMRR over 115dB

120s Real-time ECG waveforms freezing and review

Automatic measurements and interpretation tested by CSE/AHA/MIT database

Extended printing when arrhythmia is detected

Integrated advanced ECG analysis: Vector ECG



Built-in thermal recorder



8-inch LCD screen



USB support



Alphanumeric keyboard



Barcode scanner support



One-touch operation



Shortcut keys



PDF, SCP, FDA-XML, DICOM

DATA transmitting to PC via LAN/WIFI

PDF, SCP, FDA-XML and DICOM data export options

Bi-directional data communication with DMS