

## iM3 Vital Signs Monitor



Combining together a stunning ultra-slim design, large touch screen, and a rich set of features, the iM3 presents itself a next generation vital signs monitor.

- 120 hour trend review
- 1200 NIBP measurement storage
- 200 alarm storage
- Data upload via USB to EMR system
- 3 quick-action keys
- Customizable shortcuts
- Multiple working modes

**Standard Parameters:** SpO2, NIBP, PR

**Optional Parameters:** SunTech NIBP, Quick temp, Infrared Ear Temp

### Built-in Memory

Monitoring Mode (for each patient file)

<b>200</b>	<b>1200</b>	<b>120 h</b>
Alarm	NIBP Measurements	Trend Review

Spot & Round Mode

<b>16,000,000</b>
Measurements

# Innovative Design

## Adaptive Work Modes

### Round Mode

Designed to make daily use more efficient  
Easily upload data through USB or downloaded from the EMR  
Features color-coded list  
Windows interface design

### Monitoring Mode

Designed for continuous patient monitoring  
Real-time data, alarms, and history trend

### Spot Mode

Designed for spot-check applications  
Real-time vitals  
No alarm

## TEMP Module Options



Infrared Ear TEMP

Quick TEMP

Covidien Quick Temp



Portable Design



8" Touch Screen Display



360° Alarm Light



USB, OTG, LAN and DCHP



Barcode Scanner Support



Built-in WiFi

## Standard Accessories

NIBP Cuff, Adult, 27cm-35cm, reusable with NIBP tube  
Adult Finger Clip Sensor 1 meter with 2 meter extension cable  
Internal Li-Ion Battery  
Operator's Manual

## Sensor Options

SPO2 Finger Sensor, Adult, 1m, reusable  
SpO2 Wrap Sensor, Neonate, 1m, reusable  
SpO2 Silicone Soft-tip Sensor, Pediatric, 1m, reusable  
SpO2 Ear Clip Sensor, Adult/Pediatric, 1m, reusable  
SpO2 Silicone Soft-tip Sensor, Adult, 1m, reusable