

EDANUSA

iM60

Veterinary Monitor



10.4-inch high-resolution color TFT-LCD screen

Full touch screen

Customizable shortcuts for simple and intuitive operation

No-fan design reduces dust and creates a no noise environment

Defibrillation and electrosurgical interference protection

Basic parameters monitoring

Optional capnography and anesthesia monitoring available

Suitable for feline, canine and other animals

Standard Parameters: 3/5-lead ECG, HR, RESP, EDAN SpO2, EDAN NIBP, PR, 2-TEMP

Optional Parameters: Nellcor OxiMax™ SpO2, SunTech NIBP, 2-IBP, EDAN G2 CO2 (Sidestream), Respiration CO2 (Mainstream and Sidestream)

1200

NIBP
Measurement

120 h


Trend
Review

60

Alarm
Review

120 s

Frozen
Waveform





iM70

12.1" Color TFT-LCD Touch Screen



iM60


10.4" Color TFT-LCD Touch Screen




No-fan Design



Pitch Tone



Thermal Recorder



Shortcut Menu

End Tidal Carbon Dioxide Monitoring



EDAN G2 CO2 (sidestream)

Superior water trap design for accurate monitoring

iCARBTM algorithm with intelligent CO2 pseudo wave identification technology

Multiple sampling accessories as options

Respiration CO2 (mainstream/sidestream)

Plug & play module design

Dehumidification tube instead of water trap

Low sampling rate of 50ml/min


Multiple Display Modes



Standard Display



Large Font



Trend View



OxyCRG



Back to Bed



Wall Mount



Rolling Stand