



MDPRO2500

Vital Signs Monitor

Combining together a stunning ultra-slim design, large touch screen, and a rich set of features, the MDPRO2500 presents itself as a next generation monitor.

- 120 hour trend review
- 1200 NIBP measurement storage
- 200 alarm storage
- Data upload via WiFi to EMR system
- 3 quick-action keys
- Customizable shortcuts and alarms
- Multiple working modes
- 8GB internal memory

Standard Parameters: SpO2, NIBP, PR

Optional Parameters: SunTech NIBP, Covidien Oral Temp, Infrared Ear Temp, Quick Oral Temp

Built-in Memory

Monitoring Mode (for each patient file)

200	1200	120 h
Alarm	NIBP Measurements	Trend Review

Spot & Round Mode

16,000,000
Measurements

Now Integrating with **PointClickCare®**

Innovative Design

Adaptive Work Modes

Round Mode



Designed to make daily use more efficient
 Easily upload data through USB from the EMR
 Features color-coded list
 Preset and save resident names and room numbers
 Up to 1,000 names can be added at once

Monitoring Mode

Designed for continuous patient monitoring
 Real-time data, alarms, and history trend

Spot Mode

Designed for spot-check applications
 Real-time vitals
 No alarm



Portable Design



8" Touch Screen Display



360° Alarm Light



Barcode Scanner Support



USB, OTG, LAN and DCHP



Built-in WiFi

Standard Accessories

NIBP Cuff, Adult, 27cm-35cm, reusable with NIBP tube
 Adult Finger Clip Sensor 1 meter with 2 meter extension cable
 Internal Li-Ion Battery
 Operator's Manual

Sensor Options

SPO2 Finger Sensor, Adult, 1m, reusable
 SpO2 Wrap Sensor, Neonate, 1m, reusable
 SpO2 Silicone Soft-tip Sensor, Pediatric, 1m, reusable
 SpO2 Ear Clip Sensor, Adult/Pediatric, 1m, reusable
 SpO2 Silicone Soft-tip Sensor, Adult, 1m, reusable
 Disposable SpO2 Sensors, Neonate- Adult Sizes Available

TEMP Module Options



Infrared Ear TEMP
MDPRO2500_NST



Quick Oral TEMP
MDPRO2500_NST.O



Covidien Oral Temp
MDPRO2500_NST.C