

M3 Vital Signs Monitor



Outpatient



Transport



Ward



Portable Design



Infrared Ear Thermometer



M3

Vital Signs Monitor

Highlighted Features

- 5.6" High resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 8 working hours
- Nurse call

800

Alarm Review


72 h

Trend Review

30,000

NIBP Measurement Review

Communication Features



- Bi-directional communication with MFM-CMS central monitoring system
- Compatible with CMS Lite PC data management software
- HL7/AIP
- Wired LAN/WiFi Support
- USB port for flash drive or pocket router for optional WiFi
- Barcode scanner support
- Export trends as XML or PDF in the central monitoring system.

Configurations



Wall Mount



Rolling Stand



Bed Rail Clamp

Model Number	SpO2	NIBP	Quick Temp	Infrared Ear Thermometer
M3_NS	✓	✓	—	—
M3_NST.C	✓	✓	✓	—
M3_NST	✓	✓	—	✓