

EZVISION

Video Laryngoscope

A fully portable, high-performance solution for emergency care, designed with ergonomic excellence in mind. Weighing a mere .55 lbs, this cutting-edge device is equipped with an anti-fog system and a 3" HD LCD screen, ensuring clear visibility in critical situations. Additionally, it seamlessly connects to larger monitors and computers, allowing for real-time photo and video capture. With its compact size and remarkable versatility, the EZ Vision becomes an indispensable tool for emergency care providers.



.55 lbs

Light
Weight

3 HR

Continuous
Usage

3"

HD LCD
Display

Features

- 3" Wide Angle LCD Video
- Waterproof
- Portable & Lightweight
- Clear Visualization
- Real-time Capture Allows for Instant Capture of Photos and Videos
- Anti-Fog Blades with Locking Feature
- Versatile Connectivity: Can Be Effortlessly Connected to Larger Monitors and Computers
- One LCD Monitor Fits All Sizes of Blades (Disposable or Reusable)
- Disposable Handle & Blades Preventing Cross Contamination

Blade Types: Adult, Pediatric, & Neonate
(Disposable Plastic MAC 2, MAC3, MAC4 Blade)

Ergonomic Design

Comfortable grip and intuitive controls, allowing for **precise** and **efficient airway management**.



Acute



Ambulatory
Surgery Centers



EMS



Blades



Pediatric Size 2
Part Number: YLO-III B2




Medium Adult Size 3
Part Number: YLO-III B3



Large Adult Size 4
Part Number: YLO-III B4

Sample Images



Specifications

Display

Size: 3"
Resolution: 1280 x 1024px
Power: 3.7V

Laryngoscope Blade (Including Camera)

Resolution: 2 mega-pixels
Power: 3.7V
Light Source: LED 800LUX
Blade: Adult, Child, Infant

Battery

Type: Rechargeable lithium battery
Charging times: > 500 times
Charging time: <5 hours
Operating time: Lasting 300 min