

EZVISION SERIES

Video Laryngoscope

A fully portable, high-performance solution for emergency care, designed with ergonomic excellence in mind! Seamlessly connect to larger monitors and computers, allowing for real-time photo and video capture. With its compact size and remarkable versatility, the EZVision series becomes an indispensable tool for emergency care providers.



.55 lbs

Light
Weight

5 HR

Continuous
Usage

3"

HD LCD
Display

Features

- Clear visualization
- 3" Wide angle LCD video
- Portable, lightweight, & waterproof (IPX7)
- Anti-fog blades with locking feature
- Disposable blades preventing cross contamination
- Real-time capture allows for instant capture of photos and videos
- Versatile Connectivity: Can be effortlessly connected to larger monitors and computers

EZVision Blades:

Pediatric MAC 2, Medium Adult MAC 3, Large Adult MAC 4

EZVision+ Blades:

Premature Mil 00, Neonate Mil 0, Pediatric Mil 1

Ergonomic Design

Comfortable grip and intuitive controls, allowing for **precise** and **efficient airway management.**

Market Type



Acute



ASC



EMS



Primary Care



EZVision Plus
(See Page 4)

EZVision Blades



Pediatric MAC 2
Part Number: YLO-III B2




Medium Adult MAC 3
Part Number: YLO-III B3



Large Adult MAC 4
Part Number: YLO-III B4

Sample Images



Specifications

Display

Size: 3"
Resolution: 1280 x 1024px
Power: 3.7V

Laryngoscope Blade (Including Camera)

Resolution: 2 mega-pixels
Power: 3.7V
Light Source: LED 800LUX

Battery

Type: Rechargeable lithium battery
Charging times: >500 times
Charging time: <5 hours
Operating time: Lasting 5 Hours

Blades & Stands




MIL 00 (Premature)
Part Number: YLO-IB00




MIL 0 (Neonate)
Part Number: YLO-IB0



MIL 1 (Child)
Part Number: YLO-IB1



MIL 00/0 Stand (Premature/Neonate)
Part Number: YLO-IS0
INCLUDED IN PURCHASE



MIL 1 Stand (Pediatrics)
Part Number: YLO-IS1
INCLUDED IN PURCHASE

EZVision+

The EZVision+ is a video laryngoscope for premature and neonate demographics that has interchangeable stands for Miller 00, 0, & 1 sizings. Features disposable anti-fog blades for clear visibility and 3" LCD Screen.



Detachable LCD Screen
Part Number: EZVision Plus

Specifications

Display

Size: 3"
Resolution: 1600 x 1200px
Power: 3.7V

Laryngoscope Blade

Luminance: LED ≥ 150 LUX
Pixels: 2 Mega-Pixels
Power: 3.7V

Battery

Power: 3.7V
Type: Rechargeable lithium battery
Charging times: >500 times
Charging time: <4 hours
Operating time: Lasting 5 Hours