



MDPRO GUARDIAN ECG

ECG Machine



30 MIN

Long-Term
Recording

10.1 IN

Wide Capacity
Touchscreen

10 HR

Data Printout
Capability

500

ECG Data

Features

- Smart AUTO workflow
- Various report formats
- Real-time interpretation
- External memory capacity
- Save user profiles & settings
- Long-term recording up to 30 min
- 3-USB ports for spirometer connection
- HDMI output for a larger screen display
- 30 Min full disclosure with live data comparison mode
- Advanced spirometry test modes COPD, FVC, SVC, MVV
- Wireless LAN connectivity with HIPAA compliant transmission
- External USB data storage in PDF, XML, MFER, BMP, JPG, DICOM

Effective And Accurate

The MDPro Guardian is engineered to optimize workflow, enhance procedural efficiency, and improve diagnostic accuracy. Equipped with an alphanumeric keyboard and robust data security features, it supports a secure and user-friendly experience. The MDPro Guardian delivers advanced interpretation and a range of digital diagnostic modes, providing fast and accurate results to support confident clinical decision-making.

Market Type



EMS



Acute Care



Primary Care



ICU

